

Pressalit Scandinavia®, incl. UN3 univertical hinge

No. 75xxx-UN3999 with cover
No. 74000-UN3999 without cover


Design/Year

Sigvard Bernadotte, Acton Bjørn, 1975, registered trade mark

Product Description


PRESSALIT toilet seat with cover, model "Scandinavia", art. no. 75 made of colour ingrained duroplast with stainless steel hinges.

- with UN3 univertical hinge of stainless steel
- centre distance: 155 mm
- load - ring seat 240 kg
- suitable for many universal toilets
- also with institutional hinges on request (objecta seat)
- 10 year guarantee




Colour Programme xxx

With cover: 000 white, 013 bahamabeige, 050 jasmin, 052 manhattan, 273 pergamon
Without cover: 000 white



Alternative Hinges

- BN3 Universal hinge
- D43 Top mounted hinge
- B38 Top mounted hinge
- BV5 Institutional hinge, universal
- BY3 Institutional hinge, fixed



Hinge Finish

999 Stainless steel

Spare Parts

- A4003 Buffer for pushing into cover, round, grey
- A9003 Buffer with screw for seat, oval, grey
- A7220 Distance bush for seat without cover
- A7222 Clamp bush, grey
- A7250 Axle bushing
- A9107 Set of straps and screws, stainless steel



25 years anniversary in 2000 and still going strong

Product Weight	Net weight approx. 3,000 kg, incl. hinge		
Packaging	Single carton 5 single cartons in a master carton 18 master cartons (90 seats) per pallet	Dimensions: 449x381x64 mm 456x388x340 mm	Weight: 3,300 kg 17,000 kg
Material Description	<p>Toilet seat: The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises 67% urea formaldehyde resin, 28% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. The raw material from the supplier contains up to 0.03% free formaldehyde and even less following curing.</p> <p>Buffers: EVA (copolymer made of ethylene and vinyl acetate).</p> <p>Hinge links: Stainless steel (V2A) model W. no. 1.4301/M in accordance with DIN 17441/85. Composition: chrome (cr) 17-19%, nickel (ni) 8-11%, carbon <0.05%, molybdenum, manganese, silicone, sulphur, titanium.</p> <p>Mount/hinges: Stainless steel (V2A) model W. no. 1.4301/M in accordance with DIN 17441/85. Composition: chrome (cr) 17-19%, nickel (ni) 8-11%, carbon <0.05%, molybdenum, manganese, silicone, sulphur, titanium.</p>		
Suitability List	Usable on universal toilets - see www.pressalit.com		
Cleaning Instructions	<p>Use a mild soap solution to clean the seat. Please take care that the seat and hinges are not left damp. Dry off any residual water or cleaning agent with a soft cloth. Avoid contact of the seat and hinges with abrasive, corrosive or chlorine-based cleaners, as these may cause damage or lead to flash rust. Therefore, when cleaning the bowl, make sure that the seat and cover are kept in an upright position until all cleaning agent has been flushed away.</p> <p>The guarantee includes faults or defects in the material of our products within a period of 10 years. Products subject to minor technical modifications and design deviations E. & O.E.</p>		



www.pressalit.com

European Sales Agencies
Unit 10c Tanfield Lea (North) Industrial Estate
Tanfield Lea, Stanley
Co Durham, DH99UU
Tel. 01207 236622, Fax 01207 236633
info@europeansalesagencies.com
www.europeansalesagencies.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel. (+45) 8788 8788
Fax (+45) 8788 8789
pressalit@pressalit.com
www.pressalit.com

PRESSALIT®
seats