

Pressalit Dania, incl. BP5 universal hinge with quick-release
Pressalit Dania, incl. B47 universal hinge

No. 982000-BP5999 with cover
No. 982xxx-B47999 with cover
No. 981000-BU1999 without cover



Design/Year

Pressalit, 1995

Product Description

PRESSALIT toilet seat with cover, model "Dania", art. no. 982 made of colour ingrained duroplast with stainless steel hinges.

- with BP5 universal hinge with quick-release/ B47 universal hinge
- centre distance: BP5: 85-180 mm/ B47: 120-170 mm
- load - ring seat 240 kg
- with invisible hinge straps for easier cleaning
- suitable for many universal toilets
- 10 year guarantee

Colour Programme xxx


With cover: 000 white, 273 pergamon
Without cover: 000 white

Hinge Finishes


999 Stainless steel

Spare Parts


- A4013 Buffer for pushing into seat, oval, grey
- A4023 Buffer for pushing into cover, oval, grey
- A7223 End bush, barbed, long, grey
- A7224 End bush, barbed for seat without cover, short, grey
- A9193924 Set of straps with screws, chrome




BP5




BP5



B47



B47



An updated classic with invisibly mounted hinge-links - leaving the inner sides of the seat and cover with an even, smooth hygienic surface - and with its own name sort of «framed and glazed».

Pressalit Dania, incl. BP5 universal hinge with quick-release
Pressalit Dania, incl. B47 universal hinge

No. 982000-BP5999 with cover
No. 982xxx-B47999 with cover
No. 981000-BU1999 without cover

Product Weight

Net weight approx. 3,200 kg, incl. hinge

Packaging

Single carton
5 single cartons in a master carton
18 master cartons (90 seats) per pallet

Dimensions:	Weight:
449x381x64 mm	3,500 kg
456x388x340 mm	18,000 kg

Material Description

Toilet seat:

The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises 67% urea formaldehyde resin, 28% cellulose, as well as 5% minerals, pigments, lubricants and moisture content.

The raw material from the supplier contains up to 0.03% free formaldehyde and even less following curing.

Buffers:

EVA (copolymer made of ethylene and vinyl acetate).

Hinge links:

Metal version: Zinc with surface treatment. Plastic version: PA6 (polyamide).

Mount/hinges:

Stainless steel (V2A) model W. no. 1.4301/M in accordance with DIN 17441/85.

Composition: chrome (cr) 17-19%, nickel (ni) 8-11%, carbon <0.05%, molybdenum, manganese, silicone, sulphur, titanium.

Suitability List

Usable on universal toilets - see the printed seat suitability specifications.

Cleaning Instructions

Use a mild soap solution to clean the seat. Please take care that the seat and hinges are not left damp. Dry off any residual water or cleaning agent with a soft cloth. Avoid contact of the seat and hinges with abrasive, corrosive or chlorine-based cleaners, as these may cause damage or lead to flash rust. Therefore, when cleaning the bowl, make sure that the seat and cover are kept in an upright position until all cleaning agent has been flushed away.

The guarantee includes faults or defects in the material of our products within a period of 10 years.
Products subject to minor technical modifications and design deviations E. & O.E.



www.pressalit.com

European Sales Agencies
Unit 10c Tanfield Lea (North) Industrial Estate
Tanfield Lea, Stanley
Co Durham, DH99UU
Tel. 01207 236622, Fax 01207 236633
info@europeansalesagencies.com
www.europeansalesagencies.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel. (+45) 8788 8788
Fax (+45) 8788 8789
pressalit@pressalit.com
www.pressalit.com

PRESSALIT®
seats